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Operating Manual  
Manuel d'utilisation  
Instrucciones de uso  
Istruzioni per l'uso

## BeeFill™

Guttapercha Backfill-Gerät  
Guttapercha Backfill device  
Appareil d'obturation Gutta-percha  
Equipo Backfill Gutapercha  
Apparecchio per otturazione con guttaperca

## Sprachen/Languages/Langues/Idiomas/Lingua

Deutsch .....	3-13
English .....	14-24
Français .....	25-35
Español .....	36-46
Italiano .....	47-57

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Seite

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# Table of Contents

Page

Standard Symbols.....	14
Package Contents.....	15
Specifications.....	15
Certifications.....	16
Indications for Use.....	16
Contraindications .....	16
Warnings.....	17
Precautions .....	17
Adverse Reactions.....	17
Step by Step Instructions.....	18
Setting Up the Unit.....	18
Operation.....	19
Replacing a depleted cartridge .....	21
Sterilization & Maintenance.....	21
Clinical Technique.....	22
Troubleshooting .....	23
Warranty.....	24

## Standard Symbols

**Fuse Ratings**  
0.315A/110V  
0.160A/230V  
SLO-BLO

**Type BF Applied Part**

**Consult Instructions for Use**

**Class II Equipment**

**Dangerous Voltage**

**Do not re-use**

**Keep away from sunlight and heat**

**Opened packages are not replaced**

**Caution, this product contains natural rubber latex which may cause allergic reactions.**

**Cartridge contains guttapercha**

**Cannula: silver**

**Cartridge: aluminium**

**Expiry date**

**Green Dot**

**EU Authorised Representation**

## **Congratulations!**

BeeFill is engineered to provide many years of reliable service. Please read the instructions provided in this manual to receive the best and longest service from your equipment. We recommend the use of a rubber dam during all endodontic procedures.

## **Package Contents**

- Obturation control device with handpiece
- Remote power cord
- Bending tool for the cartridges
- Cleaning tool for the handpiece
- Spare cartridge nut
- Operating manual

Additional contents of the obturation kit:

- 1 Package of 10 BeeFill™ guttapercha cartridges, size 20G (90µl/0.25g)
- 1 Package of 10 BeeFill™ guttapercha cartridges, size 23G (90µl/0.25g)
- 2 Machtou-Plugger size 1-2 and 3-4

## **Specifications**

<b>Console Dimensions:</b>	11.7 cm x 18.3 cm x 8.6 cm (4.6" x 7.2" x 3.4")
<b>Handpiece:</b>	Length: 18.7 cm, Width: 2.1 cm (7.35" x 0.84")
<b>Weight:</b>	0.95 kg (2.1 lbs)
<b>Power Source:</b>	110V/60Hz, 230V/50Hz
<b>Capsule Contents:</b>	90µl / 0.25g guttapercha
<b>Capsule:</b>	Aluminium
<b>Cannula:</b>	Coin Silver
<b>Fuses:</b>	0,160A, 250 V

The appliance inlet is the main disconnect means.

### **Environmental Conditions:**

Operating Temperature:	0 to 45 C (32 to 113 F)
Storage Temperature:	-20 to 60 C (-4 to 140 F)
Relative Humidity:	5 to 95% non-condensing
Altitude:	0 to 3048 meters (0 to 10,000 feet)

**Caution:** This device has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2: 1993, EN 60601-1-2:1994, Medical Device Directive 93/42/EEC

## Certifications

Manufacturer:  
Aseptico, Inc., Woodinville, WA 98072, USA

EU Authorised Representation for CE:  
Advena Ltd., Hereford, HR4 9DQ, UK

Distributor:  
VDW GmbH, Bayerwaldstrasse 15, 81737 Munich, Germany  
Phone ++49-(0)89 627 34-0 • Fax ++49-(0)89 627 34-190  
www.vdw-dental.com

This medical device was put in circulation after Aug. 18. 2005.  
The product may not be disposed to the normal house garbage.  
The directives for the waste of electronic and electrical equipment  
must be considered.

The operating manual can be requested in other languages.

## Indications for Use

BeeFill is an obturation device with a handpiece that is used for the heating and placement of guttapercha into the root canal.

The single-use cartridges are destined for placement of the heated guttapercha into the cleaned and shaped canal.

The BeeFill device may be used only with original BeeFill cartridges!

## Contraindications

**Do not use on patients with a known sensitivity to natural rubber latex, silver or copper.**

## Warnings

Equipment not suitable for use in the presence of a flammable gas.

The guttapercha cartridges contain natural rubber latex which may cause an allergic reaction.

When replacing a cartridge during a procedure, the handpiece end cap and expended cartridge may be **hot** to the touch.

Do not use a disinfection solution for sterilization of cartridges.

Risk of fire: Do not contact the heated tip or handpiece with flammable gases or liquids.

Do not use near standing water. Dropping the unit into water may cause electric shock.

Do not immerse the handpiece or cartridges into any liquid.

## Precautions

Holding the handpiece too tightly during use, may result in a broken cannula.


Take care when replacing a cartridge during a procedure, the cartridge nut and the expended capsule may be hot to the touch.

Do not clean unit with a flammable cleaning solution.

## Adverse Reactions


Use in patients with a known sensitivity to latex, silver or copper may cause an allergic reaction. Such an allergic reaction to latex may result in swollen eyes, lips or face. It may also cause difficulty in breathing. The patient should be advised to notify you immediately if any of these symptoms occur.

## Step by Step Instructions – Setting Up the Unit

1. Unpack the console and check that the voltage selector is set to the proper voltage. Use 115V for 110-120V 60Hz voltages, and the 230V for 220-250V 50Hz voltages. To change voltage, use a flat head screwdriver to turn the voltage selector on the back of the console to the proper voltage setting. If necessary, the fuses can be replaced (see troubleshooting).
  2. Remove packaging from handpiece. Wipe down handpiece with a mild disinfecting solution – do not submerge handpiece. Set handpiece in holder on control console.
  3. Attach the power cord to the back of the console and plug into a grounded electrical outlet.
  4. Install cartridge into handpiece: Unscrew and remove the front cap from the handpiece. Insert a cartridge into the handpiece, needle facing outward. Slip cap over needle and screw cap on clockwise **lightly – do not tighten.**
- NOTE:** Take care when replacing a cartridge during a procedure, the handpiece cartridge nut and the expended capsule may be hot to the touch. **Do not remove a hot cartridge!**
5. If the cartridge doesn't fit completely into the handpiece, press the power switch on the front of the console to the 'ON' position, and press the return button on the console.  
**The plunger must be in the neutral position to accept the cartridge.**
  6. Turn the unit to the OFF position before replacing the cartridge.
  7. Use needle bending tool to place a smooth radius curve on the needle so the needle can extend to within 5 mm of the working length of the canal.
  8. Place the needle between the two raised bending posts. Gently bend the needle to the desired angle.

## Handpiece Components


## Operation

After the unit has been set up, you are ready for operation as follows:

1. Press the Power Switch on the front of the console to the 'ON' position. (The last used temperature and flow settings are retained in memory)
2. Select either Preset Button for temperature and flow rate or manually program as follows (Preset values: temperature 180° C, 60% flow rate):
  - a. Press the Temperature Increase button to the desired temperature. We suggest a temperature of 180°C. The yellow Temperature LED will illuminate indicating that the numerical LED's are displaying target temperature for approximately five seconds. While the system's temperature increases, the numerical LED window will blink until the target temperature is reached. Once the target temperature is achieved, the numerical LED window will hold steady.  
**CAUTION:** The unit must not be used at temperatures below 180°C, as this would cause damage to the motor!"
  - b. Press the Flow Rate % Increase button to the desired flow rate. The yellow Flow Rate LED will illuminate indicating that the numerical LED's are displaying flow rate for approximately five seconds. We suggest starting at 60% flow rate and then adjust from there to your preference.
  - c. To store your new settings, press either Preset button and hold until the preset LED illuminates (approximately two seconds).
  - d. Settings can be stored on both Preset buttons. This will allow you to store two distinct presets.

3. Dispense GP into patient's tooth: Engage the needle within the root canal in accordance with your preferred technique. Press any flat portion of the activation cuff to start the flow of filling material. **You will notice a short delay as the plunger in the micromotor engages and pushes guttapercha to the tip of the needle.** A small amount of guttapercha will extrude from the needle. Wipe the gutta-percha from the tip before taking the needle to the canal. Note, the indicator in the window will move as the micromotor pushes the filling material into the needle. As the material is expressed, the indicator will help you estimate how much material remains in the cartridge. Hold the handpiece lightly when placing material to allow the device to readily back out of the canal.

**NOTE:** Applying excessive downward pressure or not allowing the device to back out of the canal may result in a broken needle.

## Replacing a depleted cartridge

**NOTE: When replacing a cartridge during a procedure, the handpiece end cap and the expended cartridge are hot to the touch. Do not remove a hot cartridge!**

1. Retract the plunger in the micromotor by pressing the Return button.
2. Turn the unit off.
3. Allow the handpiece to cool for 60 seconds.
4. Unscrew and remove the handpiece cartridge nut.
5. Using the included bending tool, remove the cartridge from the handpiece.
6. Dispose of used cartridge in an appropriate biohazard container.
7. Insert a new cartridge, cannula facing out, into the handpiece end.
6. Slip cartridge nut over cannula and screw on clockwise **lightly** – do not tighten
7. Cool handpiece cartridge nut before removing any excess guttapercha.

## Sterilization & Maintenance

**Control Console** – Clean the exterior of the console by wiping with a soft cloth moistened with mild detergent or disinfecting solution. The system is designed to allow for a disposable sanitary plastic barrier to be utilized over control console and/or handpiece, which is recommended.

**Handpiece** – Wipe off the handpiece with a soft cloth moistened with mild detergent or disinfecting solution.

**DO NOT SUBMERGE** the handpiece in any fluid or spray any fluid directly on the handpiece.

**NOTE:** When wiping down the handpiece cable, gently wipe from the middle of the cable out to the handpiece and console. **Avoid gripping the cable tightly.** The handpiece cartridge nut may be steam autoclaved for 10 minutes at 135° C.

**Handpiece cartridge heater** – If guttapercha gets into the cartridge heater section of the handpiece, turn on the unit and press the return button on the console to fully retract the plunger. Allow the heater section to reach operation temperature (180° C). Turn off unit. Insert the included handpiece cleaning brush into the heating chamber. Rotate the brush several times to remove the guttapercha from the chamber.

**Cartridges** – Cartridges intended for single patient use. Prior to using the device on a patient, wipe the cannula with alcohol or a disinfectant. (Make sure that the device has been switched off and cooled down.) Store cartridges at room temperature. Do not immerse the cartridges into any liquid. Dispose of cartridge in a biohazard container. The cartridges must not be used after the expiry date!”

## Clinical Technique

1. Clean and shape the canal using your preferred technique.
2. Fit a tapered master cone into the canal and radiographically confirm it is at length.
3. Dry the canal with appropriately sized paper points. Determine final working length.
4. Trim the master cone back to the coronal foramen. Adapt plugger 3-5 mm to fixed working length. Place silicone stop.
5. Lightly coat the master cone with sealer and gently slide it to length.
6. Sever the placed master cone 5 mm from the apex. Remove excess guttapercha. Heat-soften condense and downpack the master cone to within 5 mm of canal's terminus.
7. The down pack is accomplished utilizing the vertical condensation, continuous wave or other hybrid filling techniques.
8. Before backfilling, select appropriately sized pluggers for the procedure.
9. Position the tip of the hot needle against the previously packed filling material, press the flat portion of the activating cuff on the handpiece and dispense a few millimeters of thermosoftened material into the canal.

**NOTE: Do not touch the lips, the gums or the oral mucosa with the needle or the handpiece cartridge nut, as it may get hot to the touch after longer exposure. Always use a plastic barrier for the handpiece and a rubber dam.**

10. You will notice a short delay as the plunger engages and pushes guttapercha to the tip of the needle.

**NOTE:** The indicator in the window will move as the plunger pushes the filling material into the needle. As the material is expressed, the indicator will help you estimate how much material remains in the cartridge.

11. Allow the flow of filling material to lightly back the handpiece out of the canal.

12. Use a small plugger to adapt the injected thermo-softened guttapercha into this region of the canal.
13. Sustain plugger pressure for 5 seconds on the cooling material to offset shrinkage.
14. Position the tip of the hot needle against the previously packed filling material, activate and dispense another few millimeters of material into the canal.
15. Select a larger manual plugger and step its working end circumferentially around the canal to three-dimensionally mold and tightly adapt the filling material into this region of the canal.
16. Continue the backfilling technique, in the manner described, until the canal is completely filled or stop at any point to accommodate placing a post.

## Troubleshooting

### **1. Device does not turn on:**

- a) Check that power cord is plugged into the device and the wall outlet.
- b) Unplug device and check fuse. If fuse is blown, replace with 0.160A, 230/250V Slo-Blo type fuse (0.315A for 110V). For replacement: press on the black plastic clip of the fuse, pull out and replace with a 160mA fuse.

### **2. Cannot remove handpiece cartridge nut:**

- a) Press 'RETURN' button and allow drive mechanism to back out in order to relieve pressure on cartridge nut.

### **3. Guttapercha does not flow through needle:**

- a) Verify that the needle is not kinked - avoid over bending of cartridge needle.
- b) Verify that the device has reached operational temperature for guttapercha to flow. Increase temperature if needed.

## Warranty

VDW warrants this product against defects in material or workmanship for a period of one year, from date of original invoice.

VDW's sole obligation under product warranty is (at its sole option and discretion) to repair or replace any defective component or product in part or whole. VDW shall be the sole arbiter of such action.

In the event of alleged defect under warranty, the purchaser is to notify VDW Customer Service Department promptly. Customer Service will provide instructions, usually directing that the product be returned for service. Shipment to VDW and the cost thereof is always the responsibility of the purchaser.

Accidental misuse, inappropriate installation, or failure to perform directed maintenance voids the warranty.

VDW does not assume, under this warranty, any risks or liabilities arising from the clinical use of its products, whether or not such use involves coincidental utilization of products manufactured by others.

VDW makes no warranty other than that stated above, expressed or implied.

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VDW GmbH, P.O. Box 830954, 81709 Munich/Germany