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Operating Manual

BeeFill Pack™

Guttapercha Downpack-Gerät
Guttapercha downpack device

D

GB

S

SPRACHEN/LANGUAGES

| | |
|--------------|-------|
| Deutsch..... | 3-17 |
| English..... | 18-32 |
| Svensk..... | 33-47 |

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SEITE

Symbolerklärung3

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Zertifikate5

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Warnhinweise6

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Table of Contents

Page

| | |
|--|----|
| Standard Symbols..... | 18 |
| Package Contents..... | 19 |
| Specifications..... | 19 |
| Certificates..... | 20 |
| Indications for Use..... | 20 |
| Contraindications..... | 20 |
| Warnings..... | 21 |
| Precautions..... | 21 |
| Adverse Reactions..... | 21 |
| Step-by-Step Instructions - Setting Up the Unit..... | 22 |
| Operation..... | 23 |
| Battery Operation..... | 25 |
| Replacing the fuse..... | 26 |
| Sterilization & Maintenance..... | 26 |
| Clinical Technique..... | 27 |
| Instructions..... | 28 |
| Troubleshooting..... | 30 |
| Frequently Asked Questions..... | 31 |
| Warranty..... | 32 |

Symbol description

Fuses
0.100A/115V
0.050A/230V
SLO-BLO

Dangerous Voltage

Hot Surface

Type BF Applied Part

Alternating Current

Consult Accompanying Documents

Pluggers are autoclavable

Class II Equipment

1 Package content

“ON” switch

“OFF” switch

Opened packages are not replaced

Green dot

Package Contents

- BeeFill™ Pack Device
- Remote Power Cord
- Rechargeable Battery Pack
- DFU

Additional products to be used with the device:

- BeeFill™ Pluggers:
 - Small Black ISO 40/.03
 - Medium Yellow ISO 50/.05
 - Large Blue ISO 60/.06
- BeeFill™ Thermal Response Tip

Specifications

Console Dimensions: 11.7 cm x 18.8 cm x 11.2 cm (4.6" x 7.4" x 4.4")

Handpiece: 13.2 cm L x 1.5 cm W (5.2" L x 0.6" W)

Weight: 1.30 kg (2.85 lbs)

Power Source: 230 V/50 Hz, 115 V/60 Hz

Current Rating: 230 V/50 mA, 115 V/100 mA

Fuses: 115 V: 100 mA/250 V Slow-Blow® fuse
230 V: 50 mA/250 V Slow-Blow® fuse

NOTE: The appliance inlet is the main disconnect means.

Battery Type: NiMH 9.6 V, nominal 4500 mAh

Environmental Conditions:

Operating Temperature: 10° to 40°C (50 to 104° F)

Storage Temperature: -20° to 60° C (-4 to 140° F)

Relative Humidity: 5 to 95% non-condensing

Altitude: 0 to 3,048 meters (0 to 10,000 feet)

Duty Cycle: 25%

Performance characteristics:

Temperature range: 100° to 400° C

Volume 0% to 100%

Maximum continuous heating time:

10 seconds for temperatures above 200°C.

15 seconds for temperatures at or below 200°C.

One minute in thermal response mode.

Caution: This device has been tested and found to comply with the emissions requirements of IEC 60601-1-2:2001-09. These requirements provide reasonable protection against harmful electromagnetic interference in a typical medical installation. However, high levels of radio-frequency (RF) emissions from electrical devices, such as cellular phones, may disrupt the performance of this device. To mitigate disruptive electromagnetic interference, position this device away from RF transmitters and other sources of electromagnetic energy.

RX: US LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A DENTIST

Certificates

This product is classified as a Class IIa product and bears the CE marking (CE 0459). The product complies with standards; Europe: IEC 60601-1 + A1:1991 + A2: 1995, Canada: CAN/CSA-C22.2 No. 601.1 and USA: UL 2901-1 (2nd edition).

Manufacturer:

Aseptico, Inc., P.O. Box 1548, 8333 216th Street SE, Woodinville, WA 98072, USA

EU Authorized Representative:

Advena Ltd, P.O. Hereford, HR4 9DQ, UK

Distributor:

VDW GmbH, Bayerwaldstr. 15, 81737 München, Germany

This medical device was put in circulation after Aug. 18th 2005. The product may not be disposed to the normal house garbage. The directives for the waste of electronic and electrical equipment must be considered.

The operating manual can be requested in other languages.

Indications for Use

The BeeFill™ Pack is intended to heat pluggers for warming and softening gutta percha master cones and searing off guttapercha cones. The device is also intended to heat accessory tips that subject a tooth to heat in order to determine thermal response of pulp vitality. The BeeFill™ Pack device may only be used in conjunction with the BeeFill™ pluggers and the BeeFill™ thermo response tip.

Contraindications

None known.

Warnings

Risk of fire: Do not contact the heated tip or handpiece with flammable gases or liquids.

Equipment not suitable for use in the presence of a flammable mixture with air or with oxygen or nitrous oxide.

This equipment provides ordinary protection against harmful ingress of liquids.

Do not immerse the handpiece in any fluid, or spray any fluid directly onto the handpiece.

Do not use near standing water. Dropping the unit into water may cause electric shock.

Chemclaving will not be recommended for sterilization and maintenance of the tips as this would cause corrosion.

Plugger tips become hot during use. Turn off unit and allow to cool before exchanging tips.

Battery packs may leak or explode if mishandled. Do not disassemble, incinerate, or submerge battery packs in any fluid

To avoid accidental activation, remove the battery pack prior to storing the unit for extended periods of time.

Do not ship the BeeFill™Pack unit with the accessory Plugger or Thermal Response tips attached.

NOTE: Depleted battery packs should be recycled or disposed of properly.

Precautions

A dental dam should be used with any dental procedure

Always unplug the unit before changing fuses or adjusting the voltage selector.

When replacing tips during a procedure, the plugger tips may be hot to the touch.


Do not clean with a flammable cleaning solution.

Adverse Reactions

None known.

Step by Step Instructions – Setting Up the Unit

Figure 1: View of Back Panel / Fuses


***CAUTION:** Programming port is for production use only.

1. Unpack the console and check that the voltage selector is set to the proper voltage. Use the 115 V position for 110-120 V 60 Hz voltages, and the 230 V position for 220-250 V 50 Hz voltages. To change voltage, use a flat head screwdriver to turn the voltage selector on the back of the console to the proper voltage setting. See Fig. 1
2. Remove packaging from handpiece. Clean the exterior of the handpiece with a soft cloth moistened with mild non-chlorinated detergent or disinfecting solution which is approved by the DGHM, or FDA, or bears a CE Mark. Also, do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions. Do not submerge handpiece. Set handpiece in holder on control console.
3. Attach the power cord to the back of the console and plug into a grounded electrical outlet.
4. Allow 16 hours for battery pack to fully charge before initial use.
5. Unpack the plugger and thermo tip. Sterilize pluggers and thermal response tip: Steam autoclave for 10 minutes at 132°C.
6. Install plugger into handpiece: Push plugger into handpiece tip and slowly rotate it until it keys-in with handpiece. Push plugger in fully. The device is now ready for operation.

WARNING: Do not touch hot pluggers!

Handpiece Components

Operation

After the unit has been set up, you are ready for operation as follows:

1. Turning the device on:

Press the Power switch on the bottom left-hand side of the console to the 'ON' position.

NOTE: This switch applies battery power to the unit – it does not turn AC power On and Off. The BeeFill™ Pack charges its battery whenever the unit is plugged into an AC wall outlet, regardless if the power switch is On or Off. (The last used temperature and tone level settings are retained in memory).

2. Preset Use:

Select either preset button for temperature and tone level or manually program as described in step 3, below: The selected preset temperature and tone-level values will be displayed on the numerical LED window and the preset button's LED will illuminate, indicating that its settings have been enabled. (The preset temperature is 200°C and preset tone level is set at 40%)

3. Manually program new preset temperature and tone-level settings:

a. Press the temperature up/down buttons to set the desired temperature. Press either button repeatedly to increase/decrease the temperature in 10° C increments, or hold down the button to rapidly scroll through the temperature range (100° C minimum to 400° C maximum). When a button is pressed, the temperature settings are displayed on the numerical LED window and the yellow temperature LED illuminates, indicating that the temperature setting mode has been enabled.

NOTE: The thermal response mode temperature is preset to 90° C and is not adjustable.

b. Press the volume up/down buttons to set the desired audible tone volume. Press either button repeatedly to increase/decrease the volume level in 20% increments, or hold down a button to rapidly scroll through the volume range (0% minimum to 100% maximum). When a button is pressed, the volume settings are displayed on the numerical LED window and the yellow volume LED illuminates, indicating that the volume setting mode has been enabled.

NOTE: The numerical LED window automatically defaults to displaying the selected temperature setting. If either of the volume up/down buttons is pressed, the numerical LED window will display the volume for 5 seconds and then revert to the selected temperature setting.

c. To store your new temperature and volume settings for later recall, press either of the two preset buttons and hold for approximately two seconds. The specific preset button's LED will illuminate when the settings have been successfully stored.

4. Operation:

Press the circumferential cuff on the handpiece to begin heating its tip to the selected temperature setting. If the Temperature LED is illuminated, the selected temperature will display on the numerical LED window. As the tip's temperature rises, one tone frequency (approximately 2500 Hz) is heard.


Once the tip is within 20°C of the selected temperature setting, the tone will change to a lower frequency (approximately 2000 Hz) and remain at this tone until the handpiece cuff-switch is released. When pressing the cuff, the tip will heat for a maximum time of 15 seconds if the selected temperature setting is 200°C or lower. For temperature settings above 200°C, the tip will heat for no more than 10 seconds before timing out and turning off. When in the thermal response mode, the tip will heat for up to one minute while the cuff switch is depressed. To restart heating after a timeout has been reached, release the handpiece cuff and then press it again.

Note: The temperature setting on the console has been optimized based on clinical review. Therefore, the setting does not directly correlate to the absolute temperature of the heated plugger

5. Standby Mode:

Standby mode is indicated by the decimal points of the numerical LED window illuminating in a scrolling sequence. The device will revert to standby mode after five-minutes of inactivity. If the battery charge is low (less than 10.3 volts) the unit will enter standby mode after only two minutes of inactivity. Pressing any button on the keypad will bring the unit out of standby mode and restore normal operation.

**Figure 6 –
Control Panel
Components**

Battery operation

The BeeFill™ Pack charges its battery pack whenever the unit is plugged into an AC wall outlet, regardless of the power switch position. When the battery is charging, the green LED indicator on the back cover of the device illuminates. A “standby” mode is provided to conserve battery power when the device is not in use.

NOTE: The device can't be used when drained or when no battery is installed!


Figure 2 - Aligning Battery Pack Connector


Figure 3 – Battery Pack Installed

Standby Mode:

Standby mode is indicated by the decimal points of the numerical LED Window illuminating in a scrolling sequence. The device will revert to standby mode after five-minutes of inactivity. If the battery charge is low (less than 10.3 volts) the unit will enter standby mode after only two minutes of inactivity. Pressing any button on the keypad will bring the unit out of standby mode and restore normal operation.

Low-Battery Indication:

When the battery voltage falls below 10.3 volt the device will beep three times. The numerical LED window will display “LO” for approximately five seconds and then revert to normal operation.

NOTE: When battery voltage is low, the device will enter standby mode after two minutes of inactivity. Fully charge the battery pack before performing a procedure.

Installation of new battery pack

(See Figs. 2 & 3)

- (a) Lay device on its face on a padded surface to avoid scuffing during battery installation.
- (b) Place battery pack into the cavity in the bottom of the device and align connectors as shown in Fig. 2.
- (c) Fully engage the connectors – the connector tabs will click together. Place coupled connectors into battery compartment as shown in Fig. 3.
- (d) Install the battery pack cover. Align the battery pack cover to the base of the device. By hand, align the screws with the holes and begin threading the screws into the device

chassis. Using a #2 Phillips head screwdriver, tighten down all four screws - do not overtighten.

CAUTION: Only charge using the BeeFill Pack controller.
Do not: incinerate, heat, disassemble, or short.

**Figure 4 –
Battery Pack
Cover Plate Installed**

Replacing the fuse

NOTE: The BeeFill Pack is manufactured with 250V 50mA rated fuses installed.

WARNING: Turn the power off and unplug the unit before following the steps below.

1. Remove the Fuseholder from the Power Inlet connector.


2. Replace the fuses in the Fuseholder.

Replacement Fuses:

230V: 50 mA, 250V rated slow-blow type (Fuse size: 5 x 20mm)

115V: 100 mA, 250V rated slow-blow type (Fuse size: 5 x 20mm)

3. Replace the Fuseholder.

Sterilization & Maintenance

Control Console – Clean the exterior of the console by wiping with a soft cloth moistened with mild non-chlorinated detergent or disinfecting solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives or very acidic or very alkaline solutions.

CAUTION: DO NOT SUBMERGE the device in any fluid. The system is designed to allow for a disposable sanitary plastic barrier to be utilized over control console and handpiece. Replace disposable plastic barrier between each patient use.

NOTE: When wiping down the handpiece cable, gently wipe from the middle of the cable out to the handpiece and console. Avoid gripping the cable tightly.

Handpiece – Clean the exterior of the handpiece with a soft cloth moistened with mild non-chlorinated detergent or disinfecting solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions.

CAUTION: DO NOT SUBMERGE the handpiece in any fluid or spray any fluid directly onto the handpiece.

Pluggers and Thermal Response Tip – Pluggers and tips must be cleaned, disinfected and sterilized before every use:

1. Gently use a brush for removal of debris. Clean the disinfected pluggers and tips under running water for a minimum of 1 minute. Clean with water and mild, non-abrasive, detergent or disinfection solution which is approved by the DGHM or FDA or bears a CE mark. Also do not use disinfecting solutions which contain phenol, anticorrosives, very acidic or very alkaline solutions. Clean the disinfected pluggers and tips under running water for a minimum of 1 minute to remove all chemicals. Dry thoroughly.
2. Steam autoclave plugger tips for 10 minutes at 132°C before initial use and between each patient use.
3. Observe that the pluggers and thermal response tip you are using are not damaged. A deformed or oxidized tip should be replaced. All tips will gradually lose their heating efficiency over a period of time.
4. The tips can be re-used, if they are treated with care, are not damaged or contaminated. Any further use of damaged or contaminated tips is at the user's own risk. All risk and liability is excluded in this case.

NOTE: Observe the concentration of the cleaning or disinfecting solution stated by the manufacturer!

WARNING: Do not touch hot pluggers!

Clinical Technique

Use the BeeFill™ Pack together with BeeFill™. Only use filling materials compatible with this filling technique such as Guttapercha or root canal sealer.

The actual temperature in the root canal largely depends on the amount of filling material used. Although the device is able to precisely control the heat transmission to the tip, it is unable to prevent undesired heating up of the root. We recommend testing on extracted teeth.

ATTENTION: Do not penetrate the root canal with a heated tip for more than 4 seconds for safety!

ATTENTION: Recommended temperature setting is 200° C for all downpack techniques!

NOTE: Do not touch the lips, the gums or the oral mucosa with the plugger the cannula or the handpiece cartridge nut, as it may get hot to the touch after longer exposure.

Instructions

Cone Fit & Plugger Selection

1. Prepare the canal optimally recognizing shaping facilitates 3-D cleaning and filling.
2. Select a larger size manual plugger that will work passively and effectively over a range of a few millimeters in the coronal one-third of the canal.
3. Select a medium size manual plugger that will work passively and effectively over a range of a few millimeters in the middle one-third of the canal.
4. Select a smaller size manual plugger that will work passively, effectively and deeper in the straightaway portion of the canal and to within 4-5 mm of the canal terminus.
5. Select an electric heated plugger that will passively fit through the straightaway portion of the canal and optimally to within 5 mm from the working length. Set the silicone stop at this depth to promote safety and accuracy.
6. In a fluid-filled canal, fit a non-standardized, fully tapered Guttapercha master cone that visually goes to the working length, exhibits apical tugback and is confirmed radiographically.
7. Dry the canal with appropriately sized paper points to determine final working length.
8. Trim the master cone back to the canal terminus based on the paper point drying technique.
9. Lubricate the master cone lightly with sealer and gently insert it to length.

Downpack: Continuous Wave Technique

1. Activate the electric heat plugger and sear off the master cone at the level of the orifice.
2. Select the larger manual plugger and step its working end around the circumference of the canal using short, firm strokes to move guttapercha apically, clean the canal walls and flatten the material.
3. Use this larger manual plugger and press for five seconds to compact warm guttapercha vertically and laterally into this region of the root canal system (first wave of condensation).

4. Activate the electric heated plugger and, in one deliberate and continuous motion, press the heated plugger through the thermosoftened guttapercha until the silicone stop is 2 mm from the reference point. This procedure must be limited to two to four seconds to prevent thermal injury.
5. Deactivate the plugger and continue to maintain firm apical pressure on the cooling instrument until the silicone stop reaches the reference point.
6. Maintain firm apical pressure for ten seconds to compact the mass of warm guttapercha into the apical one-third of the root canal system and to offset shrinkage during the cooling phase.
7. Activate the electric heated plugger for one second, then deactivate and remove the plugger from the root canal using a back and forth motion. This procedure separates and removes guttapercha from the coronal two-thirds of the canal without disturbing the guttapercha in the apical one-third.
8. Select the small manual plugger and step its working end around the circumference of the canal using short, firm strokes to clean the canal walls and re-condense the most coronal aspect of the guttapercha within the apical one-third.
9. Use BeeFill™ to optimally fill the rest of the canal.

Backfill

1. Position the tip of the warm cannula against the previously packed filling material for 5 seconds.
2. Press the activation cuff on the BeeFill™ handpiece and dispense a small (a few millimeters) segment of warm guttapercha into this region of the canal. Hold the handpiece lightly so it will back out of the canal during use.
3. Select the smaller manual plugger and step its working end circumferentially around the canal using short, firm strokes to clean the canal walls and flatten the dispensed warm guttapercha.
4. Use the same small manual plugger and press for five seconds to three-dimensionally compact warm guttapercha into this region of the canal and to offset shrinkage during the cooling phase.
5. Position the tip of the warm cannula against the previously packed filling material for five seconds.
6. Press the activation cuff on the BeeFill™ handpiece and dispense a few more millimeters of warm guttapercha into this region of the canal.

7. Select the medium size manual plugger and step its working end circumferentially around the canal using short, firm strokes to clean the canal walls and flatten the dispensed warm guttapercha.
8. Use the same medium size manual plugger and press for five seconds to three-dimensionally compact warm guttapercha into this region of the canal and to offset shrinkage during the cooling phase.
9. Continue the backfilling technique, in the manner described, until the canal is completely filled or stop at any point to accommodate a post to facilitate the restorative needs.

Troubleshooting

1. Device does not turn on:

- a. Verify that power cord is plugged into the device and into wall outlet.
- b. Unplug device and check fuse. If fuse is blown, replace with appropriate fuse:
115 V: 100 mA/250 V Slow-Blow® fuse or 230 V: 50 mA/250 V Slow-Blow® fuse
- c. The device has no battery installed. Always have the battery installed.

GB

2. Handpiece does not work

- a. The device has no battery installed. Always have the battery installed.
- b. Make sure the battery is charged for 16 hours. If the battery has a very low charge, plug in the power cable and wait until the battery has reached a minimum charge.
- c. Make sure the tip is sitting in the handpiece correctly.
- d. Replace tip.
- e. Replace battery.

3. The tip does not get warm

- a. Check if battery is low.
- b. Check temperatures settings.
- c. Replace tip with a new one.

Frequently Asked Questions

What size is the smallest tip?

The smallest is the small black plugger with ISO 40 Taper 03. The material of the outer surface is made out of stainless steel.

How hot can a tip get?

At maximum power, a tip can reach 400° C. You can adjust the heat down to 100° C by using the lowest power setting. The recommended working temperature for downpack is 200° C.

Isn't this high temperature dangerous?

Like other dental equipment, the unit should be used by a trained professional. You can use the high heat setting to **sear off** guttapercha in the lateral condensation technique where a cluster of cones are cut through. If you want to **remove** guttapercha, with continuous wave technique you should use a lower heat setting of 200° C.

What is the normal power setting?

The maximum heat setting should be used only when searing off guttapercha. For the vertical condensation technique and continuous wave technique, we recommend a power setting of 200°C.

How long will the tip last?

It depends on how well you take care of the tip. The life of a tip will be reduced by:

- a. running it too hot for long periods of time
- b. working with high temperatures
- c. bending or putting too much mechanical force on it.

The tips should be used only to transfer heat – not to exert manual force! Our recommendation is that you use a cold manual plugger for that purpose. All tips will gradually lose their heating efficiency over time. This is strongly depending on the user.

How long does the battery hold a charge?

The BeeFill™ Pack battery has an average operation time from 1 to 2 hours and can be recharged in 16 hours.

What is the battery life time?

The battery is without charge memory as it is a nickel metal hydride battery. However, the battery has a charge cycle life. The battery life depends on the number of charge cycles and the use of the battery. There will be no degradation in the battery life when the unit is left plugged in.

The device can't be used when drained or when no battery is installed!

Can I work with the battery only or running of the current?

Yes, you can work both while the device is plugged in charging as well as on the battery only.

Warranty

VDW warrants the BeeFill™ Pack device against defects in material or workmanship for a period of one year from date of original invoice. VDW does not warrant any other products that come along with the device, such as pluggers or thermo tips for example, as they are consumables.

VDW's sole obligation under product warranty is (at its sole option and discretion) to repair or replace any defective component or product in part or whole. VDW shall be the sole arbiter of such action.

In the event of alleged defect under warranty, the purchaser is to notify VDW Customer Complaint Department promptly. Customer Complaint will provide instructions, usually directing that the product be returned for service. Shipment to VDW and the cost thereof is always the responsibility of the purchaser.

Accidental misuse, inappropriate installation, or failure to perform directed maintenance voids the warranty.

VDW does not assume, under this warranty, any risks or liabilities arising from the clinical use of its products, whether or not such use involves coincidental utilization of products manufactured by others.

VDW makes no warranty other than that stated above, expressed or implied.

GB

Distributor:

Bayerwaldstr. 15
D-81737 Munich
Germany

Manufacturer:

Aseptico, Inc.
P.O. Box 1548
8333 216th Street SE
Woodinville, WA 98072
USA

EU Authorized
Representative:
Advena Ltd.,
Hereford, HR4 9DQ, UK

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Advena Ltd, Hereford, HR4 9DQ,