# DT LIGHT<sup>®</sup>SL


**SAFETY  
LOCK<sup>®</sup>**  
TECHNOLOGY

 **VDW<sup>®</sup>**

**Endo Easy Efficient<sup>®</sup>**

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### INDICATIONS\*

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration.

\* Sorrentino R, Monticelli F, Goracci C, Zarone F, Tay FR, Garcia-Godoy F, Ferrari M  
Effect of post-retained composite restorations on the fracture resistance of endodontically-treated teeth.  
American Journal of Dentistry, 2007, 20: 269-274

Abb 1 / Fig. 1/ Schéma n°1/ rys 1/ obr 1/ рис 1


Abb.2/ Fig.2/ Schéma n°2/ rys.2/ obr.2/ рис.2

# DT LIGHT<sup>®</sup>SL

|    |                          |    |
|----|--------------------------|----|
| de | GEBRAUCHSANWEISUNG       | 1  |
| en | INSTRUCTIONS FOR USE     | 10 |
| fr | MODE D'EMPLOI            | 18 |
| pl | INSTRUKCJA UŻYTKOWANIA   | 27 |
| cs | NÁVOD K POUŽITÍ          | 36 |
| ru | ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ | 44 |

Additional languages available upon request.

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**Caution:** For dental use only.

## DESCRIPTION

DT LIGHT®SL is a radiopaque translucent quartz fibre post in a double tapered shape. Its double taper permits a good fitting to the canal form as presented after modern treatment. This post design allows for inserting the post without unnecessary loss of tooth substance while offering optimal retention for building up the residual tooth.

The post surface is pre-conditioned with silicate ( $\text{SiO}_x$ ) and silane ( $\text{SiH}_4$ ). A protective polymer layer ensures that the silicate and the silane are not contaminated or deactivated. This polymer will be chemically integrated in the resin cement during the luting process. The Safety Lock® coating does away with pretreating the post with silane or primer and ensures a long-term bonding of the post with the composite cement.

The initial length of the post is 20 mm. The apical length of the coated surface is 15 mm. For post dimensions refer to fig.2.

## INDICATIONS FOR USE

In case of insufficient residual tooth substance (more than one missing coronal wall), the post is needed to support the coronal restoration. A ferrule effect (on at least 2 mm of coronal wall) should be possible. If not, the tooth should be treated periodontically or orthodontically in order to create a sound coronal dental structure.

## **CONTRAINDICATIONS**

Tooth with a residual circular coronal substance of less than 1.5 - 2 mm.

## **WARNING**

None known.

## **PRECAUTIONS**

- The post must be cleaned with alcohol before insertion into the canal.
- Do not touch the posts with your fingers after cleaning.
- Do not coat the post with silane.
- Shorten the post outside of the mouth.
- Remove the coloured ring from the post before inserting it in the canal.
- VDW recommends the use of a rubber dam whenever possible.

## **ADVERSE REACTIONS**

None known.

## **STEP-BY-STEP INSTRUCTIONS**

1. Preparation, cleaning and obturation of the root canal  
Prepare, clean, dry and fill the root canal according to the scientifically recognised techniques of endodontics. Avoid the use of root canal filling materials

containing eugenol. They can impair a hardening of the luting composites. VDW recommends the use of a rubber dam whenever possible.

2. Select the correct DT LIGHT®SL post according to the anatomical situation, the radio graph and the indications on fig. 1.
3. Determine the post length considering that approximately 1/2 to 1/3 of the total post length must be inserted into the root canal. Remove the root filling with the DT Universal Drill, ref. 0935 000 (rotation speed 1000-2000 rev/min.). At least 4 mm of root canal filling should remain apically. However, in the presence of severely curved canals, this calculated length cannot always be achieved. The unnecessary loss of tooth substance is to be prevented.
4. Shape the canal with the DT Finishing Drill (ref. 0935 001 or 0935 002 or 0935 003) according to the selected DT LIGHT®SL at a rotation speed between 1000-2000 rev/min. If post size 0 is used, no additional drilling should be necessary. The DT Universal Drill is the precision drill for size 0. It should be ensured that all of the gutta-percha is removed from the canal walls to facilitate an adhesive bonding – this applies in particular to oval canals.
5. Rinse the root canal (aqua dest.) to remove dentine chips and dry with paper points.
6. Clean the post with alcohol before fitting.
7. Check if the post properly fits in the canal.
8. Remove the coloured ring and shorten the coronal post to its final length outside of the mouth with a diamond disc. Never use a crimping instrument, such as wire cutters, as

the pressure can destroy the structure of the post.

9. Clean the post with alcohol. Remove excess alcohol or allow it to vaporise.

### **Adhesive bonding:**

10. DT LIGHT®SL quartz fibre posts can be bonded with all of the customary adhesive procedures and all composite cements on a BisGMA- or UDMA basis. A total etch – total bond procedure is recommended for the highest bonding values.

The following products were successfully tested with the Safety Lock® coating:

- Adhesive systems: self- or dual curing:  
Prime&Bond NT, XP Bond + SCA (Dentsply); Clearfil SE BOND (Kuraray); Syntac, Excite DSC, Adhese (Ivoclar Vivadent); Scotchbond (3M Espe); Gluma Comfort Bond (Heraeus Kulzer); ALL-BOND 2, ONE-STEP, ONE-STEP PLUS (Bisco)
- Composite cement: self- or dual curing:  
Core X-flow, SmartCem 2, Calibra, X-Flow (Dentsply); Panavia F 2.0 (Kuraray); Variolink II, Multilink Automix (Ivoclar Vivadent); LuxaCore (DMG); Rely X Unicem, RelyX ARC (3M Espe); ParaCem (Coltene Whaledent); 2Bond2 (Heraeus Kulzer); Duo Link (Bisco); G-Cem (GC)
- Core build-up: composite material

**Adhesive bonding using the example of the total etch – total bond procedure:**

11. Apply etching gel\* in the canal and on the exposed dentine in accordance with the manufacturer's instructions. Next rinse and dry gently with paper points, but leave canal wall surface moist. Avoid contact with gums.
12. Apply adhesive\* (dual- or self-curing) with a very thin brush (e.g. VDW Composibrush) in the root canal and on the coronal dentine. Remove excess with paper points and gently air dry. The surfaces should have a uniform glossy appearance. If not, repeat application and drying as described above. Light cure if needed.\*  
DT LIGHT®SL has a Safety Lock® coating, which is made of a silane and silicate layer, protected against contamination and deactivation by a polymer layer. The Safety Lock® coating does away with pretreating the post with silane or primer and ensures a long-term bonding of the post with the composite cement.
13. Mix the dual curing composite cement\* (do not use any exclusively light-curing composite cements) and apply it in the canal using a cannula, then apply composite cement on the post. Seat the post immediately. Remove excess with an appropriate instrument or a foam pellet. Afterwards light-cure\* and in the meantime stabilise the post by applying gentle pressure with the tip of the curing light fibre optic.
14. For a subsequent composite build-up of the exposed post, apply the bonding agent\* with the post (only if less than 5 mm of the quartz fibre post has been shortened) and create the build-up of the post with the composite material.\*

\* Please refer to the manufacturer's instructions for use.

A one-stage procedure is likewise possible (post insertion and post build-up in one go). It is recommended that the coronal end of the post is completely covered with the composite in order to avoid the possible penetration of liquid along the quartz fibre and thereby a hydrolytic degradation of the bonding agent.

**Retreatment:**

If a retreatment of the post should be necessary, it is recommended to use the DT Post Removal Kit (ref. 0938 001). The DT LIGHT®SL quartz fibre posts can be drilled using the help of the DT Pilot Drill and the DT Carbide Drill.\*

**STERILISING, DISINFECTING AND MAINTENANCE**

DT LIGHT®SL are single use devices. They must not be re-used to avoid the risk of contamination.

Disinfect DT LIGHT®SL with medical isopropyl alcohol. Do not use a thermo disinfectant.

**DT Universal Drill and DT Finishing Drill**

The instruments are supplied unsterile. Please disinfect and sterilise instruments before each use. Disinfect instruments with disinfecting or cleaning agents appropriate for rotary instruments and follow the manufacturer's instructions for use. Make sure instruments do not touch in the ultrasonic cleaner. Do not disinfect instruments in thermo disinfectant.

\* Observe the instructions for use of the DT Post Removal Kit.

Sterilise instruments in autoclave with the following settings\*:

- Wrap individually in disposable sterilisation pouches or tubing that comply with the regulation standard in the country of use (e.g. EN ISO 11607-1)
- Autoclave: type B complying with regulation standard in force (e.g. EN 13060)
- Sterilisation temperature: 134°C – sterilisation time: 18 min

Control instruments visually before each use: change instruments when they are damaged, corroded or dull. It is advised not to use the drill more than 12-15 times.

## **STORAGE**

Keep out of direct sunlight. Store at room temperature.

Store the sterilised components in a dry, dust-free place.

If the packaging's integrity appears to be compromised, clean and sterilise according to the protocol defined before.

\*Please refer to your autoclave manufacturer for specific sterilisation instructions if needed.

## **SAFETY AND LIABILITY**

### **Warranty**

The Manufacturer warrants this product will be free from defects in material and manufacture. The Manufacturer makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and the Manufacturer's sole obligation shall be repair or replacement of the Manufacturer product.

### **Limitation of Liability**

Except where prohibited by law, The Manufacturer will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability. Selection of the correct drill and post is the responsibility of the user.

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VDW GmbH
VDW GmbH • Bayerwaldstr. 15 • 81737 Munich • Germany