en

# FileCare® EDTA

Root Canal Conditioner

**Composition**  
FileCare® EDTA Root Canal Therapy Gel is a combination of ethylenediaminetetracetic Acid (EDTA) and carbamide peroxide in a water soluble base. EDTA cleanses calcium salts from calcified areas within the root canal to facilitate the use of reamers and files.

**Indications for Use**  
FileCare® EDTA is a device in the form of a gel for use in the chemical and mechanical cleansing of the root canal preparation during endodontic therapy. FileCare® EDTA is a combination of EDTA and carbamide peroxide specifically formulated to provide cleansing of the root canal preparation and facilitate shaping of the root canal. The FileCare® EDTA allows for cleansing action that facilitates easy removal of vital pulp tissue and necrotic pulp tissue. The FileCare® EDTA is designed to be used with endodontic irrigation with sodium hypochlorite solutions. Oxygen bubbling occurs through the release of oxygen from the carbamide peroxide. This action allows for pulp tissue, dentinal shavings and debris to be easily removed. FileCare® EDTA encourages lightening of the tooth if discoloration exists from non vitality.

The use of sodium hypochlorite with the FileCare® EDTA promotes internal bleaching of the tooth. This process is enhanced by the release of oxygen from the carbamide peroxide. FileCare® EDTA encourages lightening of the tooth if discoloration exists from non vitality. Also, the FileCare® EDTA leaves an appropriate internal canal surface for penetration of medicaments used in endodontic therapy by more completely cleaning the canal and opening dentinal tubules.

**Contraindications**  
Do not use FileCare® EDTA outside of the Root Canal System or on soft tissues. Immediately wash skin with soap and water after contact or flush eyes if contact occurs.

**Warnings:**  
1. Carbamide Peroxide and EDTA should not be applied to tissue outside of the Root Canal such as mucosa or skin. Wash skin with soap and water after contact.  
2. FileCare® EDTA canal conditioner should extrude easily; DO NOT USE EXCESSIVE FORCE. Replace original cap of FileCare® EDTA tightly after each use to avoid evaporation.  
3. Do not take internally.

**Precautions:**  
The viscosity of FileCare® EDTA is affected by storage conditions. If a more viscous or thicker consistency is desired, store in a refrigerator.

## STEP BY STEP INSTRUCTIONS

1. Open the pulp chamber in the usual manner and obtain adequate access to the root canal or canals.
  2. Remove the pulp and flush canal with sodium hypochlorite solution.
  3. Dispense the FileCare® EDTA from the tube to a dappen dish and carry to canal on an instrument such as lentulo spiral, or an endodontic instrument.
  4. The hypochlorite will react with the FileCare® EDTA and begin to effervesce or bubble.
- PLEASE TAKE NOTE OF THE BUBBLING REACTION.**
5. Enhance this reaction by mixing the FileCare® EDTA and Sodium Hypochlorite solution with a fine file.
  6. The first instrument should not be worked more than half way down the root canal.
  7. Follow the first instrument with a gentle irrigation with Sodium Hypochlorite solution.
  8. Repeat the procedure and carry the second instrument to 3/4 of the way through the canal. TAKE NOTE OF THE LUBRIFICANT OR EMULSIFYING EFFECT OF THE CONDITIONER WHEN INSTRUMENTING.
  9. Irrigate and instrument to 1 mm short of the radiographic apex.
  10. Take an X-ray to confirm measurement of the working length.
  11. The FileCare® EDTA is used in the canal with only the first two or three sizes of instruments.
  12. The canals should be further prepared using only Sodium Hypochlorite.
  13. Sodium Hypochlorite solution should be used in the canals and repeatedly instrumented into the canals until all evidence of bubbling has ceased.
  14. Dry the canal and obturate.

## Hygiene

1. Cleaning the syringe  
Replace original syringe cap prior to storing. The syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.
2. Disinfecting the syringe  
Disinfect the syringe with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

**Storage and Shelf Life**  
FileCare® EDTA should be stored at temperatures between 2° - 26°C. Refrigeration is recommended.

**Composition/Data on components:**  
Chemical characterization  
- Description: Mixture of substances listed below with nonhazardous additions.  
- Dangerous components: 60-00-4 edetic acid  
Xi; R 36 15,0 %  
124-43-6 hydrogen peroxide-urea 10,0 %  
C; R 34

**Hazard identification**  
- Hazard description:  
Xi Irritant  
- Information concerning to particular hazards to man and environment  
R 36/38 Irritating to eyes and skin.  
- Classification system  
The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

**First aid measures**  
- After inhalation  
Supply fresh air; consult doctor in case of complaints.  
- After skin contact  
Immediately wash with water and soap and rinse thoroughly.  
- After eye contact  
Rinse opened eye for several minutes under running water.  
Then consult a doctor.  
- After swallowing  
Drink plenty of water and provide fresh air. Call for a doctor immediately.

**Exposure controls and personal protection gear**  
- Additional information about design of technical facilities:  
No further data  
- Ingredients with limit values that require monitoring at the workplace:  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.  
- Additional information:  
The lists valid during the making were used as basis.  
- Personal protective equipment  
General protective and hygienic measures  
Keep away from foodstuffs, beverages and feed.  
Immediately remove all soiled and contaminated clothing  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes and skin.  
- Respiratory protection: Not required.  
- Protection of hands: Protective gloves.  
- Eye protection: Tightly sealed goggles.  
- Body protection: Protective work clothing.

**Disposal considerations**  
- Product:  
- Recommendation  
Must not be disposed together with household garbage. Do not allow product to reach sewage system.  
Uncleaned packaging:  
- Recommendation:  
Disposal must be made according to official regulations.

The instructions are available in several languages upon request.

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