




Endo
Easy
Efficient

A close-up photograph of a dentist wearing a white surgical mask and blue scrubs, holding a white EDDY irrigation device in their gloved hand. The device has a thin, curved tip. The background is blurred, showing a dental office setting.

EDDY®
Innovative
Sonic Powered
Irrigation

Step by Step
Card



Preparation

Prepare canal to a minimum ISO 25, taper .06.

Step by Step

Important:

Only activate airscaler after the EDDY® tip is inserted in canal.

Deactivate airscaler before the EDDY® tip is removed from canal.



Irrigation

1. Screw EDDY® clockwise manually into airscaler hand-piece until resistance increases slightly.



2. Fill pulp chamber with a suitable irrigation solution for endodontics, e.g. NaOCl solution, EDTA.
3. Introduce the EDDY® tip to the prepared root canal. Do not activate the airscaler.
4. In the case of an open apex, ensure the EDDY® tip is inserted only up to 1 mm short of the apex.

Airscaler **ON**



5. Once fully inserted, activate the airscaler up to 30 seconds, moving the EDDY® tip carefully up and down with vertical motions.
6. Do not exert any additional manual pressure.
7. If there is resistance, discontinue irrigation.


Airscaler **OFF**



8. Deactivate airscaler before removing the EDDY® tip from the root canal.
9. Repeat irrigation with EDDY® after every preparation step.



10. After irrigation: Remove EDDY® from airscaler hand-piece manually turning it counter-clockwise.



Obturation

Dry canal and obturate.



VDW GmbH

Bayerwaldstr. 15

81737 Munich

Germany

Tel +49 (0) 89 62734-0

Fax +49 (0) 89 62734-304

info@vdw-dental.com

