

Mtwo[®] TORQUE VALUES

Mtwo [®]	gcm	Ncm	Mtwo [®]	gcm	Ncm
10/.04	120	1.2	25/.07	200	2.0
15/.05	130	1.3	30/.06	120	1.2
20/.06	210	2.1	35/.06	100	1.0
25/.06	230	2.3	40/.06	170	1.7
30/.05	120	1.2			
35/.04	120	1.2	Retreatment instruments		
40/.04	160	1.6			
45/.04	160	1.6	R15/.05	30	0.3
50/.04	200	2.0	R25/.05	120	1.2
60/.04	300	3.0			
All sizes:			280 rpm		

Files can be operated between 250-350 rpm. The VDW torque values are recommended for safety and efficiency.

100 gcm are equivalent to 0.981 Ncm.

