

High-tech apex locator for precise length determination

Sets new standards with respect to user-friendliness and design




- Accurate length determination
- Foldable, pocket-sized design
- Unique 3-D style colour touch screen
- Smart user interface
- "Plug-and-start" and "plug-and-check" features provide visual cues for enhanced ease of use
- Rechargeable battery



Precise length determination


- Latest multi-frequency apex locator technology
- Automatic calibration
- Separate RAYPEX Apical Zoom shows the enlarged section between apical constriction and apical foramen
- Visualisation of the apical constriction

Coloured RAYPEX Apical Zoom



User friendly interface

- Unique colour touch screen for smart operation
- Intuitive menu for setting options




- DEMO mode for clear information to patient
- Built-in CHECK mode for smooth functional test of device and cables: Automatic "plug-and-check" with visual cues



Customisation for individual preference

- Choice between 2 sounds and 2 background display colours



- Easy setting of optional Dr's Choice apical line



Range of products

- REF 1113 000 000 RAYPEX®6 Kit
 1 RAYPEX®6 apex locator
 1 measuring cable
 2 lip clips
 2 touch probes
 1 file clip
 1 charger
 1 tester for functional check

Accessories

- REF 0141 000 514 5 lip clips
 REF 0141 000 515 2 file clips
 REF 0141 000 518 2 touch probes
 REF 1119 000 506 1 measuring cable
 REF 1119 000 501 1 charger
 REF 1119 000 507 1 rechargeable battery