

Cutting edge technology: The new smart apex locator VDW.CONNECT Locate[®] is soon available

Munich, March 12th, 2019 | Easy-to-control and highly accurate – VDW expands its Endo-System for root canal treatment with its innovative VDW.CONNECT Locate[®]. Combined with the smart and cordless VDW.CONNECT Drive[®] endo motor and the VDW.CONNECT[®] App, the new apex locator offers a variety of features that support the workflow of clinicians. This includes simultaneous monitoring of file progression while shaping, which allows live feedback during the preparation of the root canal.

VDW.CONNECT Locate[®] is a smart apex locator, which provides clinicians an easy and reliable way to detect the file position for their root canal treatments. The device is a smart apex locator with a state-of-the-art user interface, which connects via Bluetooth[®] with the VDW.CONNECT[®] App. This way file progression can be visualized simultaneously on an iPad screen. Additional acoustic signals and LEDs directly on the apex locator support the safe detection of the file. The device is portable and during the treatment placed on the tray.

VDW.CONNECT Locate[®] can be optionally used with the VDW.CONNECT App[®]. Only this way, clinicians benefit from additional features: The integrated shaping target function enables the clinician to customize until which point to shape the root canal. The apical reverse function, if enabled, automatically rotates the file back in the opposite direction of the cutting direction when the previously defined reference point is reached.

This maximum flexibility, great ergonomics of the system and numerous functions enable the clinician to perform a safe treatment.

Discover more about VDW.CONNECT Locate[®] at IDS 2019 in hall 10.1, B050 or at www.vdw-dental.com.

Contact

Christiane Rohn
Marketing Communication Manager

VDW GmbH
Bayerwaldstr. 15
81737 Munich
Germany

T +49 (0)89 627 34 - 114
F +49 (0)89 627 34 - 195
christiane.rohn@vdw-dental.com

ABOUT VDW:

VDW GmbH based in Munich, Germany is one of the most well-known manufacturers working in the dental field of endodontics. For 150 years, VDW has been a pioneer in shaping the evolution of root canal treatment significantly. VDW focuses on offering the dentist a holistic solution covering the entire endodontic treatment process including preparation, irrigation, obturation and post-endodontic care as well as service and training. For more information about the company, the VDW brand and products, please visit www.vdw-dental.com.

Pictures:

Pic. 1: VDW.CONNECT Locate[®] - a smart, portable apex locator for efficient endo procedures.



Pic. 2: The VDW.CONNECT[®] family integrates the new apex locator VDW.CONNECT Locate[®], the smart and cordless handpiece VDW.CONNECT Drive[®] and the VDW.CONNECT[®] App.

