

## Identification of the product

Product name: Silver points  
Use of the article: root canal post

### Manufacturer/supplier identification:

Company: VDW GmbH  
Bayerwaldstr. 15, D-81737 München, Germany Tel: ++49 89 62734-0  
e-mail: info@vdw-dental.com  
Contact for information: Quality Management  
Telephone: ++49 89 62734-0 • Telefax: ++49 89 62734-304

## Composition / information on ingredients

Root canal post of a Silver - Nickel special alloy.

	contains
Silver	80-100%
Nickel	< 0,2%

## Physical and chemical properties

Diameter:	0,55 und 0,75 mm
Expansion	2%
Tensile strength:	about 395 N/mm <sup>2</sup>

## Other information

The information contained herein is based on the present state of our knowledge. It characterizes the technical information of the product.