



RECIPROC[®] direct
The first reciprocating
contra-angle

Make your
dental unit an
endo unit

Unique and efficient: Experience the first reciprocating contra-angle

- RECIPROC® direct contra-angle connects directly to your existing dental unit 
- Simply click it on your electric or air-driven motor* and get access to RECIPROC® one file endo
- Profit from innovative engineering that drives RECIPROC® instruments safely and efficiently



Plug and play: Prepare for your treatment in three easy steps



Attach the contra-angle to the motor so that it snaps into place



Make sure RECIPROC® direct is securely attached to the motor



Set your motor* on your dental unit to the appropriate speed

Air-driven motor: max. speed

Electric motor: 20.000 rpm (max. 25.000 rpm)

*ISO-E coupling according to ISO 3964 is necessary

Innovative and affordable: Your starter kit with RECIPROC® direct

Best value for money:

Contra-angle + 24 RECIPROC® *blue* instruments

RECIPROC® direct Starter Kit with RECIPROC® *blue*

1 x RECIPROC® direct contra-angle

Blister à 6 instruments 25 mm

● 2 x Blister of RECIPROC® *blue* instruments R25

● 1 x Blister of RECIPROC® *blue* instruments R40

● 1 x Blister of RECIPROC® *blue* instruments R50



Order number: V04 1095 025 000



For further information please contact your preferred distributor or visit www.vdw-dental.com

VDW GmbH
Bayerwaldstr. 15
81737 Munich
Germany
Tel +49 (0) 89 62734-0
Fax +49 (0) 89 62734-304
info@vdw-dental.com

