























# EDDY<sup>®</sup>

## Airscaler Kompatibilität

Airscaler		Kupplung	
3H <sup>®</sup>	 Sonic S SS-NP	NSK oder Phatlus Kupplung	
NSK	 Ti-MAX S970, S970L (Optic)	NSK Kupplung	
	 Ti-MAX S970KL (Optic)	Kavo <sup>®</sup> LED Kupplung	
	 Ti-MAX S970SL (Optic)	Sirona <sup>®</sup> LED Kupplung	
W&H	 ZA-55, ZA-55 L (Optic)	W&H Roto Quick RQ 34	
	 ZA-55M, ZA-55 LM (Optic)	Kavo <sup>®</sup> Multiflex <sup>®</sup>	
	 ZA-55 LS (Optic)	Sirona <sup>®</sup> LED Kupplung	
Kavo <sup>®</sup>	 SONICflex 2003	Kavo <sup>®</sup> Multiflex <sup>®</sup> oder Sirona <sup>®</sup> LED	
	 SONICflex 2003 L		
	 SONICflex 2000N		
Sirona <sup>®</sup>	 SIROAIR L	Sirona <sup>®</sup> R/F Quick	
Micron Corp	 TA-200 KL	Multiflex	

Für weitere Informationen kontaktieren Sie Ihren bevorzugten Händler oder besuchen Sie uns auf [www.vdw-dental.com](http://www.vdw-dental.com).