

VDW.FLO™ Endo Organizer

Multifunctional drawer insert,
equipment trolley and table
tray in a single unit

M00VFLO000000



Ideal overview

The VDW.FLO™ Endo Organizer brings order and creates a system for your materials. Organize loose instruments, irrigation, obturation and post endodontic materials to keep track of your inventory.

Practical helper

Use the VDW.FLO™ Endo Organizer as a tray and have everything you need at hand during treatment.

Customizable

You can assign each space individually with the enclosed labels or use the preprinted stickers to adapt the VDW.FLO™ Endo Organizer specifically to your workflow.



Product description



390 x 300 x 38 mm

- Plastic insert for drawers, equipment trolleys and table tray
 - suitable for standard drawer, equipment trolley and table sizes
 - 16 small, medium and large compartments matched to standard packaging for a flexible division
 - robust, long-lasting, easy to clean plastic
 - stackable
- 14 silicone feet for a non-slip stand
- 88 writable, wipeable and removable plastic stickers